


OEM Solutions for Anesthesiology and Respiratory Support



In biomedical signal processing, gas monitoring and respiratory support since 1989



Content

OEM Solutions for Anesthesiology and Respiratory Support

Depth of Anesthesia and Sedation Module	3
Multigas Module	9
Mainstream CO ₂ Sensor QuRe®	13
Sidestream CO ₂ Sensor	17
Indirect Calorimetry Module	22

By developing advanced medical solutions, we are helping doctors save patients' lives

About

Triton Electronic Systems develops and supplies customized solutions for medical device industry since 1989. We focused on gas monitoring solutions and biomedical signal processing.

Our strategy is to provide ready-to-use and plug&play solutions for our customers, to support your developments and to reduce your leadtime for launching new functions to the market.

Our main values

- Fast and easy integration.
- Highest protection (IP, EMC, safety).
- Regulatory & documentary support.
- Original design and private labeling.
- Cutting edge patented technologies.

Inspired & expired gas monitoring (CO₂, O₂)

Respiratory monitoring

Volatile anesthetics monitoring

Vital signs monitoring

Cardiac output monitoring

Sedation monitoring

Gross floor area

more than 9000 sqm

Number of employees

260 persons

R&D engineers

60 persons


Global sales network

more than 40 countries

Depth of Anesthesia and Sedation Module



Complete solution for anesthesiologists:
Monitoring of Sedation Level




Depth of Anesthesia and Sedation Module is designed to provide long and continuous monitoring of the Brain Activity Index (AI).

Application: anesthesiology, including perioperative period, resuscitation, intensive care, procedural sedation.

This is the solution for the daily routine depth of anesthesia monitoring, a standard monitoring tool in a medical institution, thereby increasing the patient's safety and quality of patient care.

Depth of Anesthesia and Sedation Module can be connected to the monitoring host device and transfers parameters of the Brain Activity Index, Signal Quantity Index, Electromyographic Component, Suppression Rate and additional states and statuses.




Identified Parameters

AI	Brain Activity Index	Indicates the level of consciousness depression by analyzing EEG, taking into account information on typical signs of anesthetics' impact on patients
SR	EEG Signal Suppression Rate	Reflects the relative duration of EEG suppression segments in the current time interval
SQI	EEG Signal Quality Index	Reflects noise influence on EEG signal
EMG	Electromyographic Component Level	Indicates the level of electrical activity of facial muscles
Z1, Z2, Z3	Electrode impedance	Demonstrates the quality of electrodes application and electrodes' electric contact with the patients' skin

Interpreting AI (Brain Activity Index) Data*

AI value	Clinical stages of anesthesia	Clinical signs
90–100	Awake	
80–90	Anesthesia stage I — light sedation	Incomplete awakening, patient opens eyes and maintains visual contact in response to a voice for 10 seconds or less
60–80	Anesthesia stage II — sedation	Patient moves and opens eyes in reaction to voice but does not fix the eyes — no visual contact or no response to voice but eye movements and eye opening after a physical stimulation persists
40–60	Anesthesia stage III — surgical state	No response to voice or physical irritants
30–40	Anesthesia stage IV — deep anesthesia, BS (burst-suppression) patterns emerge	
20–30	Anesthesia stage V — deeper anesthesia compared to stage IV, length of suppression episodes may reach 10 seconds	
0–10	Anesthesia stage VII — extremely deep anesthesia, suppression episodes constitute 75% and more of the whole signal duration	

* According to a generally accepted classification of anesthesia stages.



Advantages of Anesthesia Depth Monitoring

Potential Effects of Inadequate Sedation*

With continuously raising requirements to ensuring patient's safety, physicians have to provide more careful control of using anesthetics, hypnotic drugs or sedatives. According to the statistics, more than 69% of patients demonstrate inadequate sedation — either insufficient or too deep. This can cause adverse effects both during the surgery and at the post-operative stage.

Potential Effects of Inadequate Sedation*

Insufficient sedation	Excessive sedation
Excitation	Depressed breathing, hypotonia, depressed gastrointestinal tract motility
Sleep violations	Prolonged depression of consciousness
Myocardial ischemia	Prolonged ventilation duration
Unsynchronized ventilation	Prolonged stay at ICU and clinic in general
Self-extubation	Increased healthcare costs

Posttraumatic distress and depression

* Mehta S. Sedation Strategies in the Critically Ill // Yearbook of Intensive Care and Emergency Medicine, 2005.

Using MGA minimizes the adverse effects of inadequate sedation, ensuring optimal and predictable sedation level and patient's quicker recovery from anesthesia.

Preventing Anesthesia Awareness

Anesthesia awareness is post-operative recollections of the event happening during general anesthesia, caused by misalignment between the need for anesthetic and its delivery.

The following patients are in the risk group for anesthesia awareness:

- taking opiates or alcohol;
- using neuromuscular relaxants;
- suffering from respiratory problems;
- with previous cases of accidental awakening during the surgery;
- with a co-pathology;
- elderly.

Anesthesia awareness cannot be measured directly. Traditional clinical signs like motions, tachycardia, hypertension, pupillary reaction and lacrimation are supposed to be unreliable predictors of anesthesia awareness but they must be monitored in every patient and considered substantially.

Individual Selection of Sedative Doses


Selecting the optimal drug dosage with MGA is based on EEG analysis and displaying the AI (Brain Activity Index), taking into account individual body features and the clinical situation.

This approach ensures:

- maximum safety and efficiency of the delivered anesthetic support;
- reducing the risk of drugs' adverse effect on the body;
- saving expensive drugs.

Algorithm of AI (Brain Activity Index) Calculation

Anesthesia depth assessment is based on a comprehensive electroencephalogram (EEG) analysis using unique algorithms developed by Triton Electronic Systems engineers. A simplified algorithm for AI (Brain Activity Index) calculation is performed below.



EEG and EMG signals are registered from the electrodes applied on the frontotemporal area of the patient's head.

The registered signal is subjected to digital filtration: motion artifacts, power main disturbances and noise from electrosurgical equipment, other bioelectrical signals, etc. are removed.

The algorithm of EEG analysis includes statistical information on typical signs of various groups of anesthetics' impact on the patient's EEG. During the analysis, the level of compliance is established between the registered EEG signal and each type of conscience depression.

As a result of data analysis, the following indicator values are obtained:

AI (Brain Activity Index);
 SR — EEG Signal Suppression Rate, taking into account the total duration of segments with low-voltage activity (suppression segments) over the last minute. Displayed as a percentage. SR > 0 is usual at AI < 50.

Signal quality

To get accurate data on anesthesia depth monitoring:

- assess the signal quality continuously;
- provide control of electrodes' impedance;
- prevent impedance values from increasing;
- minimize artifacts and other noise.

For this purpose, the following technical solutions are implemented in MGA Module:

SQI (Signal Quality Index) is continuously monitored. It takes into account the values of EEG cable electrode' impedance, noise level from artifacts, high-frequency noise and power main disturbances within EEG, etc.

If SQI = 0, displaying the values of AI (Brain Activity Index), SR and EMG Component Level is blocked. A message on the most significant cause for SQI dropping is transferred.

The level of EEG signal noise is measured continuously.


Electrode impedance is an important parameter to indicate a quality of contact between skin and electrodes. Low impedance is an essential condition to get high quality signal.

Technical Specification

Patient age groups	Adults and children over 10 years old
Anesthetics	Works with both inhaled and intravenous anesthetics
Displayed parameters	Brain Activity Index EEG Signal Suppression Rate EEG Signal Quality Index Electromyographic Component Level Electrodes impedance (Z1–Z3)
The recommended types of electrodes	31.1245.21, 24 mm, Covidien LLC, USA; F9079/RU3236-100, FIAB SpA, Italy; White Sensor 40713, Ambu A/S, Denmark; G210C/F-150S, Nihon Kohden Corporation, Japan
Dimensions & Weight	115x65x25 mm, 0.2 kg
Power	Voltage: 5.0 V±5% DC. Power consumption: 2 W
Environment	0–40°C, RH 40–80% (at air temperature 25°C), 600–800 mmHg
Integration	UART interface, ODU connector
Standards	Developed in accordance with IEC 60601-1, IEC 60601-1-2, IEC 60601-2-26, IEC 60601-2-40

OEM Delivery Kit

Depth of Anesthesia and Sedation Module MGA	TESM.943129.007-02	1 pcs.
Electrodes	100 pcs.	



Multigas Module




Complete solution for anesthesiologists:
Anesthesia Gas Monitoring



Appearance

Multigas is anesthesia gas analyzer intended for continuous sidestream measurement and monitoring of gas concentration in the patient's airways. Specially designed for implementation into medical devices, including anesthesia machines, patient monitors and ventilators.



Multigas Module

Measured Parameters

Concentration of anesthetic gas agents (Ax): Isoflurane (Iso), Sevoflurane (Sev), Desflurane (Des) on inspiration (FiAx) and expiration (EtAx), Halotane (Hal) — option.

Concentration of CO₂ on inspiration (FiCO₂) and at the end of exhalation (EtCO₂).

Respiration rate (RR).

Concentration of O₂ (option) on inspiration (FiO₂) and at the end of exhalation (EtO₂).

Multigas benefits

- Easy integration.
- Maintenance-free.
- Standard accessories.
- No routine calibration.



Agent Measuring Ranges

Gas	Measurement range, vol. %	Accuracy
Isoflurane (Iso)	0–5	±(0.2% + 15% of gas level)
Sevoflurane (Sev)	0–7	±(0.2% + 15% of gas level)
Desflurane (Des)	0–17	±(0.2% + 15% of gas level)
Halothane (Hal) (option)	0–5	±0.2
CO ₂	0–15	±(0.43% + 8% of gas level)
O ₂ (option)	0–100	±2

OEM Delivery Kit

Multigas Module	TESM.943129.001	1 pcs.
Water trap	60-1300-00 DRYLINE	1 pcs.
Sampling line	010-700 Flexicare	1 pcs.
Connector	15M-22M/15 REF2713 or 15M-22M/15 010-638 Flexicare	1 pcs.
Interface cable	TESM.704021-03	1 pcs.
Developer kit includes: RS-232-USB converter, USB cable, CD with development software and data protocol description, integration manual	On request ET555494	1 pcs.

Technical Specification

Technology	Non-dispersive infrared (NDIR)
Measured gases	Anesthetic agents: Iso, Sev, Des, Hal (option), CO ₂ , O ₂ (option)
Measured parameters	Inspired and expired gas concentrations of all gases; respiration rate; ambient pressure
Patient age groups	Adult, pediatric
Self-diagnosis	Availability
Warm-up time	ISO accuracy within 45 s. Full accuracy within 10 min, CO ₂ channel 10 s
Response time	2.5 s
Resolution	0.1
Gas sampling rate	50–250 ml/min
Respiration rate range	0–160 breath per minute (BPM)
Respiration rate accuracy	±1 breath
Apnea detection range	10–60 s, default 20 s
Calibration	No user calibration required Supports automatic atmospheric pressure compensation
Dimensions & Weight	External module: 150x95x60 mm. Weight: 0.5 kg
Power	Voltage: 5.0 V±5% Power consumption: 2.5 W, not more than 5 W at the warm-up
Environment / Protection	Water and splash resistance: IP22 Operation: 10–35°C, RH 10–90% Storage: 5–40°C, RH < 80% at 25°C, 390–900 mmHg Transportation: (–50)°C–50°C
Integration	Interface: RS-232 Connector: 5-pin ODU Data output: FiCO ₂ , FiO ₂ , FiAx, EtCO ₂ , EtO ₂ , EtAx, respiration rate
Standards	Developed in accordance with the requirements: EN 60601-1, IEC 60601-1-2, EN ISO 80601-2-55, MDD 93/42/EEC, RoHS

Mainstream CO₂ Sensor QuRe[®]



Mainstream CO₂
monitoring



QuRe® Sensor is intended to measure FiCO₂ and EtCO₂ concentration in the mainstream breathing. Designed to use with medical devices for intensive care, ventilation, respiratory support and patient monitoring. The sensor is applicable for all intubated patients from adults to neonates.




QuRe® Sensor Benefits

- The most accurate capnogram.
- Compact and lightweight.
- Long-life time reusable airway adapters, up to 100 sterilization cycles.
- Suitable for high-frequency ventilation, up to 200 BPM.
- No user calibration required.
- “Everything on board”: the measurement channel and analysis CPU are inside the sensor. This saves space in your own medical devices.
- Easy integration to any monitoring or ventilation system, full technical support by the Manufacturer.
- Private labeling solutions (PLM) available.
- Fastest response time.


OEM Delivery Kit

Mainstream CO ₂ Sensor	TESM.506001	1 pcs.
Reusable airway adapter adult / pediatric pediatric / neonate	TESM.706020 TESM.706021	1 pcs.
Evaluation kit includes: RS-232-USB converter, USB cable, CD with development software and data protocol description, integration manual	On request	1 pcs.

The highest precision of capnogram waveform for true clinical diagnostics.



Mainstream CO₂ Sensor



QuRe[®] data rate is 100 Hz and instant accurate response of the sensor to the changing CO₂ concentration.

The sensor contains an innovative ultrafast semiconductor emitter with a high modulation frequency of the signal.

Due to the precise capnography waveform, the clinicians can analyze each breathing cycle, detect respiratory disorders and provide optimal ventilation mode to a patient.

Technical Specification

Sensor type	CO ₂ mainstream sensor
Operation principle	Non-dispersive infrared (NDIR)
Initialization time	20 s at an ambient temperature of 25°C Full specification within 2 min
Response time	~10 ms
CO ₂ measurement range	0–20% (0–150 mmHg)
CO ₂ accuracy	0.2 vol.%
CO ₂ resolution	0.1 vol.%
Respiration rate range	0–200 breath per minute (BPM)
Respiration rate accuracy	±1 breath
Apnea detection range	10–60 s, default 20 s
Calibration	No user calibration required
Airway adapters	Reusable adult / pediatric. Reusable pediatric / neonate Material: airway adapter — polycarbonate, optical windows — sapphire glass. Dead space: <5 ml (adult / pediatric), <1 ml (pediatric / neonate) Sterilization: ethylene oxide / autoclaving up to 100 times
Dimensions & Weight	Sensor: 38x35x23 mm. Weight: 28 g. Cable length: 3 m
Power	Voltage: 5.0 V±5%. Power consumption: 1.3 W, not more than 3 W at the warm-up
Environment / Protection	Water and splash resistance: IP44 (sensor) Operation: 10–35°C, RH 10–90%, 390–900 mmHg Storage: 5–40°C, RH < 80% at 25°C Transportation: (–50)°C–50°C
Integration	Interface: RS-232. Connector: Lemo Redel / ODU Data output: FiCO ₂ , EtCO ₂ , respiration rate. Gas and pressure compensation: available, supplied by host

Sidestream CO₂ Sensor



Easy
integration



Intended Use

Sidestream CO₂ Sensor is intended for continuous noninvasive monitoring of fraction of inspired CO₂ (FiCO₂) and end tidal CO₂ (EtCO₂), respiration rate (RR) and apnea. Specially designed for implementation into medical devices for anesthesiology, intensive care, as well as patient monitors.

Scope Anesthesiology and intensive care departments of professional medical facilities, transportation within professional medical facilities.

CO₂ monitoring is recommended by various international associations as a routine procedure for patients during anesthesia and in intensive care departments.

Benefits

- Easy integration.
- Maintenance-free.
- Standard accessories.
- No routine calibration.

Clinical Application

Assessing of the spontaneous breathing adequacy. Using the capnography the level of spontaneous breathing during recovery after anesthesia can be assessed.

Weaning from mechanical ventilation.

Controlling of breathing system hermetic seal. A gas leakage is always possible during anesthesia. The leakage can be detected by EtCO₂ monitoring, which value is gradually increasing due to hypoventilation.

Examination of circulatory arrest and resuscitation procedures effectiveness.

Capnometry is optimal method for monitoring of cardiopulmonary resuscitation (CPR) effectiveness.

Control during the trachea intubation.

Examination of ventilation-perfusion ratio mismatch.

Any cause that reduces the lungs perfusion and/or increases the respiratory

dead space can lead to P_{Et}CO₂ decreasing and ΔP(a-Et)CO₂ increasing.

Monitoring of hypermetabolic conditions (malignant hyperthermia, thyroid crisis, sepsis, etc.).

References

1. ASA Standards for Basic Anesthetic Monitoring, Standards and Practice.
2. European Resuscitation Council Guidelines for Resuscitation.
3. American Heart Association (AHA) Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care.
4. Association of Anaesthetists of Great Britain and Ireland. Recommendations for standards of monitoring during anaesthesia and recovery.

Appearance



Sidestream CO₂ Sensor

OEM Delivery Kit

Sidestream CO ₂ Sensor	TESN.943129.008-01	1 pcs.
Connector	15M-22M/15 REF2713 or 15M-22M/15 010-638 Flexicare	*
Monitoring line	TESN.943129.008-01	1 pcs.
Water trap	60-1300-00 DRYLINE	1 pcs.
Interface cable	TESM.704021-03	1 pcs.

* The items and quantity of the accessories and documents are specified at the order.

Technical Specification

Operation principle	Non-dispersive infrared spectrophotometry (NDIR)
Initialization time	10 s
Time of setting the operating mode	120 s (at ambient temperature 25°C)
Gas sampling rate accuracy	50–250 ml/min ±10 ml/min in absolute terms or ±10% in relative one (the biggest from the values)
CO ₂ concentration measurement range	0–20 vol.% (resolution 0.1) 0–150 mmHg (resolution 0.1)
CO ₂ measurement accuracy	$\pm(0.2 \text{ vol.\%} + 0.02 \cdot K_{\text{meas}})$ or $\pm(1.5 \text{ mmHg} + 0.02 \cdot K_{\text{meas}})$
Measurements drift	$\pm(0.2 \text{ vol.\%} + 0.02 \cdot K_{\text{meas}})$ or $\pm(1.5 \text{ mmHg} + 0.02 \cdot K_{\text{meas}})$
Response time	~3 s
Rise time	0.2 s
Influence of gas impurities and vapors	$\pm(0.43 \text{ vol.\%} + 0.08 \cdot K_{\text{meas}})$
Respiratory rate measurement range Respiratory rate measurement accuracy	3–160 breath per minute (BPM) ±2 breath per minute (BPM)
Power	Voltage: 5.0 V±5% Capacity: 1.5 W, maximum 4 W during warming
Weight	0.5 kg
Dimensions (without cable), width x height x depth	58x92x146 mm
Connector	Lemo Redel / ODU
Interface	RS-232

Operation Conditions

Ambient temperature	10–35°C
Relative humidity	10–90% (at the temperature 25°C)
Ambient pressure	600...800 mmHg

Transportation Conditions

Ambient temperature	(–50)°C–50°C
Relative humidity	<80% (at the temperature 25°C)

Storage Conditions

Ambient temperature	5–40°C
Relative humidity	<80% (at the temperature 25°C)

Standards

The device meets the safety requirements of EN 60601-1, EN ISO 80601-2-55.

The device is developed to meet the requirements of EN 60601-1-2 concerning electromagnetic compatibility (EMC).

Indirect Calorimetry Module



Metabolic needs
evaluation



Intended Use

The device is intended for continuous non-invasive monitoring of carbon dioxide (CO₂) concentration and resting energy expenditure (REE) of ventilated patients.

Scope anesthesiology, resuscitation and intensive care departments of professional medical facilities, as well as at transportation within professional medical facilities. Applicable as a part of anesthesia and respiratory equipment and medical devices for patient monitoring.

Appearance

LED indicator

Water trap

Sampling line




Interface cable

OEM Delivery Kit

Calorimeter module	TESN.650002	1 pcs.
Interface cable	TESN.704029	1 pcs.
Water trap	DRYLINE II	1 pcs.
Sampling line	DRYLINE	1 pcs.
Paracube® Sprint oxygen sensor	00502755	Option

Evaluation of Patient's Metabolic Needs



The peculiarity of patients in intensive care and resuscitation units is metabolic instability caused by the severity of the condition, artificial lung ventilation, sedation, analgesia and extracorporeal detoxification methods. Therefore, metabolic monitoring for such patients is of great importance.

The method of indirect calorimetry is considered the “gold standard” of metabolic monitoring. In addition to actual resting energy expenditure

(REE), this method calculates the respiratory quotient (RQ) — the ratio of carbon dioxide release rate to oxygen consumption rate and assess the contribution of each macronutrient to the total metabolism.

The built-in metabolic module is convenient and easy to use because it requires minimal user effort.

The principle of the metabolic needs evaluation is based on measuring the volume of carbon dioxide released, the

volume of oxygen absorbed and the subsequent calculation of energy costs using the Weir equation.

Experience has shown that the individualized program of nutritional support for 3–4 days of treatment in ICU using the metabolic module significantly reduces:

frequency of nosocomial infections; consumption of antibacterial drugs; duration of artificial ventilation.

Metabolic monitoring is used in programmes of early and resuscitation rehabilitation of patients. Its use makes it possible to shorten the time of rehabilitation and minimize complications after suffering strokes, spinal cord injuries, brain injuries, muscular dystrophies, etc.

Deficiency of Calories in Critical States Can Cause

- postoperative wound suppuration, failure of anastomoses;
- dysfunction of the respiratory musculature and diaphragm;
- hospital-acquired infections (tracheobronchitis, VAP, etc.);
- high consumption of antibiotics;
- greater consumption of blood components (FFP, albumin);
- pressure sores, anemia;
- prolonged bed rest in ICU and inpatient department.

Excess Calories in Critical States Lead To

- hyperglycemia;
- growth of CO₂ production;
- desynchronization with the ventilator;
- hyperthermia;
- aggravation of ALI / ARDS;
- fatty hepatosis.

Operation Principle

The Indirect Calorimetry Module provides continuous measurement of CO₂ partial pressure in patient's airway by infrared spectrophotometry. The method consists in measuring the absorption of infrared radiation with a certain wavelength by CO₂ molecules which is calculated on the basis of measured amount of light transmitted through the gas to the device. The oxygen concentration measures by O₂ paramagnetic sensor.

The real time flow data at the patient side (Y-piece) is transmitted to the module by the host device (ventilator). The sampling probe from the patient is fed into the measurement cuvette of the module by the

integrated pump via sampling line, goes through the CO₂ cuvette (capnograph) and O₂ sensor. The module synchronizes the O₂, CO₂ and flow data in realtime by the patented algorithm to precise breath by breath measurements of minute volume of CO₂ elimination (VCO₂), and minute volume of O₂ consumption (VO₂) by the patient. The RQ (respiratory quotient) and REE (resting energy expenditure) are calculated by Weir equation: REE (kcal/day) = (3.941 x VO₂ + 1.106 x VCO₂) x N.

The essence of this method is to calculate a respiratory quotient (RQ) and a ratio of eliminated CO₂ to consumed O₂ per minute. The calculated indirect calorimetry parameters are then averaged and transferred via communication protocol to host-device (ventilator) for displaying.

Technical Specification

Operation principle	Non-dispersive infrared spectrophotometry (NDIR)		
External interface	RS-232		
Initialization time	Maximum 10 s		
Gas sampling rate Admissible absolute deviation	50–250 ml/min ±10 ml/min in absolute terms or ±10% in relative one (the biggest from the values)		
CO ₂ concentration measurement range (CO ₂ partial pressure) Admissible absolute deviation: CO ₂ concentration CO ₂ partial pressure	0–20 vol.% (resolution 0.01) 0–150 mmHg (resolution 0.01) ±(0.2 vol.% + 0.02 · K _{meas}) ±(1.5 mmHg + 0.02 · K _{meas})		
O ₂ concentration measurement range Admissible absolute deviation	0–100 vol.% (resolution 0.1) ±(2.5 vol.% + 0.025 · K _{meas})		
Calculated REE parameters: O ₂ consumption (VO ₂) CO ₂ elimination (VCO ₂) Resting energy expenditure (REE) Respiratory quotient (RQ)	10–1000 ml/min 10–1000 ml/min 72–7200 kcal/day 0.5–2		
Measurement drift	±(0.2 vol.% + 0.02 · K _{meas}) or ±(1.5 mmHg + 0.02 · K _{meas})		
Gas mixture parameter measurement accuracy	±(0.43 vol.% + 0.08 · K _{meas})		
Response time	Maximum 3 s (at sampling rate 250 ml/min)		
Rise time	Maximum 0.2 s (at sampling rate 250 ml/min)		
Respiratory rate measurement range	3–160 breath per minute (BPM)		
Respiratory rate measurement accuracy	±2 breath per minute (BPM)		
Dimensions	90x98x160 mm		
Voltage Current	5.0 V ± 5% maximum 0.8 A	Device weight	0.7 kg

Operation Conditions

Ambient air temperature	10–35°C
Relative humidity	10–80% (at the temperature 25°C)
Atmospheric pressure	600...800 mmHg

Transportation Conditions

Ambient air temperature	(–50)°C–50°C
Relative humidity	< 80% (at the temperature 25°C)

Storage Conditions

Ambient air temperature	5–40°C
Relative humidity	< 80% (at the temperature 25°C)

Export Geography



Over 40 countries around the world

Great Britain

Indonesia

Syria

Germany

Malaysia

Egypt

Portugal

China

Turkey

Spain

Taiwan

Peru

Denmark

Sri Lanka

Uzbekistan

Ireland

Bangladesh

Kyrgyzstan

Hungary

Korea

Croatia

Czech Republic

Vietnam

Belarus

Poland

Philippines

Slovakia

Thailand

Moldova

Singapore

India

South Africa



We continuously improve the technological principles and implement new profitable solutions based on market demands



In biomedical signal processing, gas monitoring and respiratory support since 1989

Bajova str. 33
Ekaterinburg, 620133
Russian Federation

Quality management system certified as meeting the requirements of EN ISO 13485

January 2024

